

Quantum Imaging and Therapeutic Associates

Pelvis (Perirectal Fistula) – Toshiba Aquilion One / 64

Indications	Perirectal fistula
Position / Landmark	Supine / iliac crest
Scan Type	Helical
Start / End Locations	Noncontrast: Iliac crests / lesser trochanters Rectal contrast: Iliac crests / lesser trochanters
FOV	28 cm (or as small as feasible to target rectum)
kV / mA	120 kV / Sure Exposure
Slice Thickness/Interval	3 mm x 3 mm
IV Contrast Volume / Rate	None
Scan Delay / Sure Start	N/A
Oral Contrast	None
MPR	Coronal and sagittal rectal contrast 3 mm x 3 mm
Images to PACS	Scouts, axial noncontrast soft tissue, axial rectal contrast soft tissue, MPR
Comments / Special	Rectal contrast for second scan. NO oral contrast.