

Quantum Imaging and Therapeutic Associates

Skull – Toshiba Aquilion One / 64

Indications	Skull injury, craniosynostosis
Position / Landmark	Supine / outer canthus of eye
Scan Type	Helical
Start / End Locations	Chin / 1 cm above vertex
FOV	25 cm (or as small as feasible)
kV / mA	120 kV / Sure Exposure
Slice Thickness/Interval	5 mm x 5 mm
IV Contrast Volume / Rate	None
Scan Delay / Sure Start	N/A
Oral Contrast	None
MPR	3D volume render of skull including all sutures
Images to PACS	Scouts, axial brain, axial bone, MPR
Comments / Special	N/A