

Quantum Imaging and Therapeutic Associates

MRV Abdomen – 1.5 Tesla

Indications	Stricture or thrombosis of IVC or iliac veins.
Coil	Flex, torso
Position	Supine
FOV	32 cm
Comments / Special	

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
2D TOF						
Oblique axial recon (aorta)					8/2	
MIP 12 projections @ 15 degrees						