

Quantum Imaging and Therapeutic Associates

Abdomen and Pelvis Screen – 1.5 Tesla

Indications	Abdominal pain, abdominal mass (unknown origin), lymphadenopathy
Coil	Flex, torso
Position	Supine
FOV	32 cm
Comments / Special	

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Cor T2 SSFSE						
Ax dual echo T1 in/out phase						
Ax T2 SSFSE						
Ax T1 FS SPGR abdomen/pelvis						
Ax T1 FS SPGR abdomen (20, 60, 120 s)						
Ax T1 FS SPGR post pelvis						