

Quantum Imaging and Therapeutic Associates

Adrenal Mass – 1.5 Tesla

Indications	Adrenal mass, suspected pheochromocytoma
Coil	Flex, torso
Position	Supine
FOV	32 cm
Comments / Special	Ideally, noncontrast portion of exam can be checked by radiologist to see if the exam is diagnostic or if contrast must be administered.

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Cor T1 SPGR						
Ax dual echo T1 in/out phase						
Ax T2 SSFSE						
Ax T2 FS SSFSE through ao bifurcation						
Ax T1 FS SPGR						
Ax T1 FS SPGR post gad (20s, 60s, 120s, 12m delay)						