

Quantum Imaging and Therapeutic Associates

Liver Mass (Eovist) – 1.5 Tesla

Indications	Liver mass.
Coil	Flex, torso
Position	Supine
FOV	32 cm
Comments / Special	

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax dual echo T1 in/out phase						
Ax T1 FS SPGR						
Ax T1 FS SPGR post gad (20s, 60s, 120s)						
Cor T2 SSFSE						
Ax T2 SSFSE						
Ax 2D SSFP FS						
Ax T1 FS SPGR post gad (20m delay)						
Cor T1 FS SPGR post gad (20m delay)						