

# Quantum Imaging and Therapeutic Associates

## Pelvis – 1.5 Tesla

Indications	Pelvic bone pain, tumor/mets
Coil	Body
Position	Supine, feet neutral/parallel
FOV	30-40 cm
Comments / Special	Optional Sag T2 FS small FOV through symphysis pubis only if for sportsman hernia

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax T1	400-800	Min.			4/1	256 x 256
Ax T2 FS	>2000	50-60			4/1	256 x 192
Cor T1	400-800	Min.			4/1	256 x 256
Cor STIR	>2000	20-40	150		4/1	256 x 192
*Sag T2 FS	>2000	50-60			3/1	256 x 256