

Quantum Imaging and Therapeutic Associates

Renal Mass – 1.5 Tesla

Indications	Kidney mass, hematuria, post-ablation follow up.
Coil	Flex, torso
Position	Supine
FOV	32 cm
Comments / Special	

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax dual echo T1 in/out phase						
Ax T2 SSFSE						
Cor T2 SSFSE						
Ax T1 SPGR						
Cor T1 SPGR						
Cor T1 FS SPGR post gad (20s, 60s, 120s)						
Ax T1 FS SPGR post gad (180s)						
Subtraction images for all postcontrast sequences						