

Quantum Imaging and Therapeutic Associates

Sacrum/SI Joints – 1.5 Tesla

Indications	Fracture, sacroiliitis, infection
Coil	Body
Position	Supine
FOV	20 cm
Comments / Special	Coronals are along axis of sacrum *Optional post gad

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax T1	400-800	Min.			4/1	256 x 256
Ax T2 FS	>2000	50-60			4/1	256 x 192
Cor T1	400-800	Min.			4/1	256 x 256
Cor STIR	>2000	20-40	150		4/1	256 x 192
Sag T2 FS	>2000	50-60			4/1	256 X 192
*Ax/Cor T1 FS post gad	400-800	Min.			4/1	256 x 256