

Quantum Imaging and Therapeutic Associates

Wrist – 1.5 Tesla

Indications	TFCC tear, tendon pathology, occult fracture, mass
Coil	Wrist, flex M
Position	Palm down
FOV	8-10 cm
Comments / Special	Skin marker at site of pain, palpable lesion, scar

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax T1	500-600	10-15			3/0.3	256 x 256
Ax PD FS	>2000	35-45			3/0.3	256 x 256
Cor T1	500-600	10-15			3/0.3	256 x 256
Cor PD FS	>2000	35-45			3/0.3	256 x 256
Cor 3D T1FFE	30-40	Min. (in phase)		30	1/0	256 x 256
Sag STIR	>2500	15-25	120		3/0.3	256 x 256