

Quantum Imaging and Therapeutic Associates

Breast (Implants) – 1.5 Tesla

| | |
|--------------------|---|
| Indications | Breast implant integrity |
| Coil | Breast |
| Position | Prone |
| FOV | Variable to fit breasts and anterior chest wall |
| Comments / Special | |

| <u>Sequence</u> | <u>TR</u> | <u>TE</u> | <u>TI</u> | <u>Flip</u> | <u>ST/sp</u> | <u>Matrix</u> |
|-----------------|-----------|-----------|-----------|-------------|--------------|---------------|
| Sag T2 | 5000 | 120 | | | 4/0 | |
| Ax T2 H2O sat | 5000 | 120 | | | 4/0 | |
| Ax IR | 3000-4000 | 15-25 | | | 4/0 | |
| Ax IR H2O sat | 3000-4000 | 15-25 | | | 4/0 | |
| Ax T1 | 500-600 | <10 | | | 4/0 | |